

Lecture 5

X-ray Photoemission Spectroscopy (XPS)

5. Photoemission Spectroscopy (XPS)


- 5.1 Principles
- 5.2 Interpretation
- 5.3 Instrumentation
- 5.4 XPS vs UV Photoelectron Spectroscopy (UPS)
- 5.5 Auger Electron Spectroscopy (AES)
- 5.6 Quantitative Analysis

References:


- 1) Zangwill; Chapter 2, pp. 20-24 and 4
- 2) Kolasinski, Chapter 2.6
- 3) Vickerman, Chapter 2 and 3
- 4) Woodruff, and Delchar, Chapter 3
- 5) Briggs, Seah, *Practical Surface Analysis*. 1991; Vol. 1.
- 6) Luth, Chapter 6
- 7) <http://www.phy.cuhk.edu.hk/course/surfacesci/index2.html>
- 8) <http://www.chem.qmul.ac.uk/surfaces/scc/>
- 9) <http://www.cem.msu.edu/~cem924sg/>

10/3/2010

Lecture 5



Electron beam interactions with the sample




10/3/2010

Lecture 5

2

9.1 Electron Scattering



Electron diffraction and microscopy:

Elastic backscattered e^- , ~ few % at 100eV

II: Interaction with plasmons

III: Inelastic electrons, Auger

IV: Secondary e^- 's, <50eV





Fig. V.2. Qualitative large-scale overview of the energy distribution of electrons emitted from a surface which is irradiated by an electron beam of primary energy E_0 .

“Universal curve” for electrons



Short inelastic mean free path for electrons means that elastic scattering of electrons is very surface sensitive


Lecture 5

3

Electron beam-solid interactions

Secondary electrons (SEs): are produced by the interactions between energetic e^- 's and weakly bonded valence e^- 's of the sample

Auger electron: incident e^- kicks out an inner shell e^- ; a vacant e^- state is formed; this inner shell vacant state is then filled by another e^- from a higher shell, and simultaneously the energy is transferred to another e^- that leaves the sample



Characteristic X-rays: emitted when a hole is created in the inner shell of an atom in the specimen due to inelastic e^- scattering, as it can recombine with an outer shell e^- (EDX)

Backscattered electrons (BSEs): are primary e^- 's leaving the specimen after a few large angle elastic scattering events

Cathodoluminescence (CL): light emission arising from the recombination of e^-h pairs induced by excitation of e^- 's in the valence band during inelastic scattering in a semiconducting sample

10/3/2010

Lecture 5

4

Electron Spectroscopy for Surface Analysis

| Spectroscopy | Particles involved | Incident Energy | What you learn |
|------------------------------|--|-----------------|---------------------------------|
| XPS X-ray Photoemission | X-ray in e out | 1-4 keV | Chemical state, composition |
| UPS UV Photoemission | UV photon e out | 5-500 eV | Valence band |
| AES Auger Electron | e in, e out; radiationless process, filling of core hole | 1-5 keV | Composition, depth profiling |
| IPS Inverse Photoemission | e in photon out | 8-20eV | Unoccupied states |
| EELS Electron Energy Loss | e in e out | 1-5 eV | Vibrations |

10/3/2010

Lecture 5




5

5.1 Photoemission Spectroscopy: Principles

Electrons absorb X-ray photon and are ejected from atom

Energy balance:

$$\begin{aligned} \text{Photon energy} & - \text{Kinetic Energy} = \text{Binding Energy} \\ h\nu & - KE = BE \end{aligned}$$



- Spectrum – Kinetic energy distribution of photoemitted electrons
- Different orbitals give different peaks in spectrum
- Peak intensities depend on photoionization cross section (largest for C 1s)
- Extra peak: Auger emission


10/3/2010

Lecture 5

6

Photoemission Spectroscopy: Basics

Electrons from the sample surface:




$$I(d) = K \int_0^d \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx$$

Fraction of signal from various depth in term of λ

| Depth | Equation | Fraction of signal ($\theta=0$) |
|------------|---|-----------------------------------|
| λ | $\frac{I(\lambda)}{I(\infty)} = \frac{\int_0^\lambda \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}{\int_0^\infty \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}$ | 0.63 |
| 2λ | $\frac{I(2\lambda)}{I(\infty)} = \frac{\int_0^{2\lambda} \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}{\int_0^\infty \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}$ | 0.86 |
| 3λ | $\frac{I(3\lambda)}{I(\infty)} = \frac{\int_0^{3\lambda} \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}{\int_0^\infty \exp\left(\frac{-x}{\lambda \cos \theta}\right) dx}$ | 0.95 |

1. C. J. Powell, A. Jablonski, S. Tanuma, et al. *J. Electron Spectrosc. Relat. Phenom*, 68, P. 605 (1994).
2. D. F. Mitchell, K. B. Clark, W. N. Lennard, et al. , *Surf. Interface Anal.* 21, P. 44 (1994).

5.2 Typical XPS (ESCA) spectrum



8

X-ray and spectroscopic notations

| | Quantum numbers | | | X-ray suffix | X-ray level | Spectroscopic Level |
|--|-----------------|-----|-----|--------------|----------------|---------------------|
| | n | l | j | | | |
| Principle quantum number: $n = 1, 2, 3, \dots$ | 1 | 0 | 1/2 | 1 | K | 1s _{1/2} |
| Orbital quantum number: $l = 0, 1, 2, \dots, (n-1)$ | 2 | 0 | 1/2 | 1 | L ₁ | 2s _{1/2} |
| Spin quantum number: $s = \pm 1/2$ | 2 | 1 | 1/2 | 2 | L ₂ | 2p _{1/2} |
| Total angular momentum: $j = l + s = 1/2, 3/2, 5/2$ | 2 | 1 | 3/2 | 3 | L ₃ | 2p _{3/2} |
| | 3 | 0 | 1/2 | 1 | M ₁ | 3s _{1/2} |
| | 3 | 1 | 1/2 | 2 | M ₂ | 3p _{1/2} |
| | 3 | 1 | 3/2 | 3 | M ₃ | 3p _{3/2} |
| | 3 | 2 | 3/2 | 4 | M ₄ | 3d _{3/2} |
| | 3 | 2 | 5/2 | 5 | M ₅ | 3d _{5/2} |
| Spin-orbit split doublets | Etc. | | | Etc. | Etc. | Etc. |


| Sub-shell | j values | Area ratio |
|-----------|------------|------------|
| s | 1/2 | - |
| p | 1/2; 3/2 | 1: 2 |
| d | 3/2; 5/2 | 2: 3 |
| 10/3/2010 | 5/2; 7/2 | 3: 4 |

Lecture 5

9

Binding energy reference in XPS

Energy level diagram for an electrically conductive sample grounded to the spectrometer



10

Typical spectral features

Associate binding energies with orbital energies, BUT USE CAUTION!


Energy conservation: $E_i(N) + h\nu = E_f(N-1) + KE$
 $\Rightarrow h\nu - KE = E_f(N-1, k) - E_i(N) = E_B$

Binding energy is more properly associated with ionization energy.

In HF approach, Koopmans' Theorem: $E_B = E_k$ (orbital energy of k^{th} level)
 Formally correct within HF. Wrong when correlation effects are included.

ALSO: Photoexcitation is rapid event
 \Rightarrow sudden approximation

Gives rise to chemical shifts and plasmon peaks



10/3/2010

Lecture 5

11

Qualitative results

- A: Identify element
 - B: Chemical shifts of core levels:
- Consider core levels of the same element in different chemical states:
 $\Delta E_B = E_B(2) - E_B(1) = E_K(2) - E_K(1)$

Often correct to associate ΔE_B with change in local electrostatic potential due to change in electron density associated with chemical bonding ("initial state effects").

Peak Width:



ΔE measured as FWHM; if assume Gaussian:

$$\Delta E = (\Delta E^2_{\text{natural}} + \Delta E^2_{\text{photons}} + \Delta E^2_{\text{analyzer}})^{1/2}$$

Lifetime width: $\Gamma = \hbar/\tau$ for $\Gamma \sim 1$ eV, $\tau \sim 10^{-15}$ s


<http://www.lasurface.com/database/liaisonxps.php>

10/3/2010



Lecture 5

Chemical Shifts



- Carbon 1s chemical shifts in ethyl trifluoroacetate
- The four carbon lines correspond to the four atoms within the molecule

10/3/2010

Lecture 5


13

Peak Identification: Core level binding energies

XPS Photoelectron Binding Energies* versus Atomic Number

Where the same BE is quoted for different spin-orbit split subshells, the individual subshells are either not resolved or have a complex structure. Spin-orbit coupling: $2p_1=2p_{1/2}$, $2p_3=2p_{3/2}$ etc.


| Z | 1s | 2s | 2p ₁ | 2p ₃ | 3s | 3p ₁ | 3p ₃ | 3d ₃ | 3d ₅ | 4s | 4p ₁ | 4p ₃ | 4d ₃ | 4d ₅ | 4f ₅ | 4f ₇ | 5s | 5p ₁ | 5p ₃ | |
|----|------|-----|-----------------|-----------------|----|-----------------|-----------------|-----------------|-----------------|----|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|----|-----------------|-----------------|--|
| 1 | 14 | | | | | | | | | | | | | | | | | | | |
| 2 | 21 | | | | | | | | | | | | | | | | | | | |
| 3 | 56 | | | | | | | | | | | | | | | | | | | |
| 4 | 113 | | | | | | | | | | | | | | | | | | | |
| 5 | 191 | | | | | | | | | | | | | | | | | | | |
| 6 | 287 | | | | | | | | | | | | | | | | | | | |
| 7 | 402 | | | | | | | | | | | | | | | | | | | |
| 8 | 531 | 23 | | | | | | | | | | | | | | | | | | |
| 9 | 686 | 30 | | | | | | | | | | | | | | | | | | |
| 10 | 863 | 41 | 14 | 14 | | | | | | | | | | | | | | | | |
| 11 | 1072 | 64 | 31 | 31 | | | | | | | | | | | | | | | | |
| 12 | | 90 | 51 | 51 | | | | | | | | | | | | | | | | |
| 13 | | 119 | 74 | 74 | | | | | | | | | | | | | | | | |
| 14 | | 153 | 103 | 102 | | | | | | | | | | | | | | | | |
| 15 | | 191 | 134 | 133 | 14 | | | | | | | | | | | | | | | |
| 16 | | 229 | 166 | 165 | 17 | | | | | | | | | | | | | | | |
| 17 | | 270 | 201 | 199 | 17 | | | | | | | | | | | | | | | |
| 18 | | 319 | 243 | 241 | 22 | | | | | | | | | | | | | | | |
| 19 | | 378 | 296 | 293 | 33 | 17 | 17 | | | | | | | | | | | | | |
| 20 | | 439 | 330 | 347 | 44 | 25 | 25 | | | | | | | | | | | | | |
| 21 | | 501 | 407 | 402 | 53 | 31 | 31 | | | | | | | | | | | | | |
| 22 | | 565 | 464 | 458 | 62 | 37 | 37 | | | | | | | | | | | | | |
| 23 | | 630 | 523 | 515 | 69 | 40 | 40 | | | | | | | | | | | | | |
| 24 | | 698 | 586 | 577 | 77 | 46 | 45 | | | | | | | | | | | | | |
| 25 | | 770 | 652 | 641 | 83 | 49 | 48 | | | | | | | | | | | | | |



<http://www.lasurface.com/database/elementxps.php>

Quantification of XPS

Primary assumption for quantitative analysis: ionization probability (photoemission cross section) of a core level is nearly independent of valence state for a given element
 \Rightarrow intensity \propto number of atoms in detection volume



$$I_A = \sigma_A(\hbar\omega)D(E_A) \int_{\gamma=0}^{\pi} \int_{\phi=0}^{2\pi} L_A(\gamma) \int_{x,y} J_0(x,y)T(x,y,\gamma,\phi,E_A) \int_x N_A(x,y,z) \exp\left(\frac{-z}{\lambda \cos \theta}\right) dx dy dz d\gamma d\phi$$

where:

- σ_A = photoionization cross section
- $D(E_A)$ = detection efficiency of spectrometer at E_A
- $L_A(\gamma)$ = angular asymmetry of photoemission intensity
- γ = angle between incident X-rays and detector
- $J_0(x,y)$ = flux of primary photons into surface at point (x,y)
- T = analyzer transmission
- ϕ = azimuthal angle
- $N_A(x,y,z)$ = density of A atoms at (x,y,z)
- λ_M = electron attenuation length of e's with energy E_A in matrix M
- θ = detection angle (between sample normal and spectrometer)




Figure 5.13 Geometry of the XPS analysis configuration

10/3/2010

Lecture 5

15

Quantitative analysis

For small entrance aperture (fixed ϕ, γ) and uniform illuminated sample:


$$I_A = \sigma_A(\hbar\omega)D(E_A)L_A(\gamma_i)J_0N_A\lambda_M(E_A)\cos\theta_iG(E_A)$$

Angles γ_i and θ_i are fixed by the sample geometry and $G(E_A) = \iint_{x,y} T(x,y,E_A) dx dy$

$G(E_A)$ = product of area analyzed and analyzer transmission function

$D(E_A)$ = const for spectrometers operating at fixed pass energy

σ_A : well described by Scofield
 Calculation of cross-section




L 5.14 Calculated values of the cross-section $\sigma_A(h\nu)$ for Al K α radiation in terms of the C 1s cross-section. (After Scofield¹⁸)

10/3/2010

5.3 Photoemission Spectroscopy: Instrumentation

X-ray source



10/3/2010

| X-ray lines | | |
|-------------------|------------|-----------|
| Line | Energy, eV | Width, eV |
| Ti L _α | 395.3 | 3.0 |
| Cu L _α | 929.7 | 3.8 |
| Mg K _α | 1253.6 | 0.7 |
| Al K _α | 1486.6 | 0.85 |
| Ti K _α | 4510.0 | 2.0 |



How to choose the material for a soft X-ray source:

1. the line width must not limit the energy resolution;
2. the characteristic X-ray energy must be high enough to eject core electrons for an unambiguous analysis;
3. the photoionization cross section of e in different core levels varies with the wavelength of the X-ray, a suitable characteristic X-ray wavelength is crucial to obtain a strong enough photoelectron signal for analysis.

Instrumentation

Essential components:

- **Sample:** usually 1 cm²
 - **X-ray source:** Al 1486.6 eV; Mg 1256.6 eV
 - **Electron Energy Analyzer:** 100 mm radius concentric hemispherical analyzer (CHA); vary voltages to vary pass energy.
 - **Detector:** electron multiplier (channeltron)
 - **Electronics, Computer**
 - **Note:** All in **ultrahigh vacuum** (<10⁻⁸ Torr) (<10⁻¹¹ atm)
 - **State-of-the-art small spot**
- ESCA: 5 μm spot size
Sputtering gun for profiling



10/3/2010


18

Electron Energy Analyzers

(a) Concentric Hemispherical Analyzer (CHA) and (b) (Double Pass) Cylindrical Mirror Analyser (CMA)

Advantages of CHA:
higher resolution than CMA, convenient geometry

Disadvantages: small solid angle



Advantages of CMA: very large solid angle

Disadvantages: inconvenient geometry (b.)


10/3/2010

Lecture 5

19

New electron energy analyzers with lens system

$E_p/\delta E \sim 1400$
 $E_p = 2\text{eV}, < 2\text{meV}$
 $\delta\Theta \sim 10\text{ mrad}$
 $\delta S \sim 9\ \mu\text{m}$



10/3/2010

Lecture 5

20

Surface Science Western- XPS

<http://www.uwo.ca/ssw/services/xps.html>
<http://xpsfitting.blogspot.com/>
<http://www.casaxps.com/>
<http://www.lasurface.com/database/elementxps.php>



Kratos Axis Ultra (left)
Axis Nova (right)


Contact:
Mark Biesinger
biesingr@uwo.ca

10/3/2010

Lecture 5

21

5.4 Comparison XPS and UPS



XPS: photon energy $h\nu=200-4000$ eV to probe core-levels (to identify elements and their chemical states).

UPS: photon energy $h\nu=10-45$ eV to probe filled electron states in **valence band** or adsorbed molecules on metal.

Angle resolved UPS can be used to map band structure (to be discussed later)

UPS source of irradiation: He discharging lamp (two strong lines at 21.2 eV and 42.4 eV, termed He I and He II) with narrow line width and high flux

Synchrotron radiation source
continuously variable photon energy, can be made very narrow, very intense, now widely available, require a monochromator

Introduction to Photoemission Spectroscopy in solids, by F. Boscherini
http://amscampus.cib.unibo.it/archive/00002071/01/photoemission_spectroscopy.pdf

Studies with UV Photoemission


- The electronic structure of solids -detailed angle resolved studies permit the complete band structure to be mapped out in k -space
- The adsorption of molecules on solids-by comparison of the molecular orbitals of the adsorbed species with those of both the isolated molecule and with calculations.
- The distinction between UPS and XPS is becoming less and less well defined due to the important role now played by synchrotron radiation.

10/3/2010


Lecture 5

23


5.5 Auger Electron Spectroscopy (AES)



Electron collision



Auger electron emission





- Steps in Auger deexcitation
- Note: The energy of the Auger electrons do not depend on the energy of the projectile electron in (a)!

10

Lecture 5

24

Auger spectrum of Cu(001) and CuP



Use dN/dE (derivative mode) \Rightarrow
Why?

Sections of Auger electron spectra, showing Cu (M2,3VV) and P (L2,3VV) transitions, for a low temperature PH₃ overlayer phase at 140K and (b) for a P c (6×8) structure obtained by annealing the surface of (a) to T_x > 450K. Both spectra have been normalized to give the same Cu (60 eV) feature peak height.

10/3/2010

Lecture 5

25

Applications of AES


- A means of monitoring surface cleanliness of samples
- High sensitivity (typically ca. 1% monolayer) for all elements except H and He.
- Quantitative compositional analysis of the surface region of specimens, by comparison with standard samples of known composition.
- The basic technique has also been adapted for use in :
 - Auger Depth Profiling : providing quantitative compositional information as a function of depth below the surface (through sputtering)
 - Scanning Auger Microscopy (SAM) : providing spatially-resolved compositional information on heterogeneous samples (by scanning the electron beam over the sample)

10/3/2010


Lecture 5

26



5.6 Quantitative analysis



The big picture: Elemental Identification



Zoom in: Chemical shifts, plasmons



Quantitative Analysis

- Estimate chemical concentration, chemical state, spatial distribution of surface species
- Simplest approximation is that sample is in single phase




FIG. 1. Idealized surface structures: (a) plane homogeneous surface; (b) a surface with lateral inhomogeneities consisting of several different surface phases; (c) a surface with depth inhomogeneities (the circles and the crosses represent different types of atoms); (d) a surface phase consisting of a submonolayer of foreign atoms on an otherwise homogeneous surface; and (e) an interface between two homogeneous bulk phases (Ref. 1).

10/3/2010

J. Vac. Sci. Technol. A, Vol. 4, No. 3, May/June 1986